

NEW LABIAL BOW

TESTED OVER 5 YEARS ON MORE THAN 6000 ORTHODONTIC APPLIANCES




Aesthetic

Adjustable

Guaranteed*

Designed to provide both aesthetics and durability, the new ClearBow is a 3 prong plier adjustable labial bow offering 2.75 mm of full labial contact. ClearBow features a patented design and is fabricated from food-grade materials made in the United States and manufactured in Canada. The BPA-free plastic strap is guaranteed not to break under normal use, and is not damaged by boiling water or stained by coffee, cola or red wine.



Superior aesthetics

Encourages patient compliance

Adjustable

Plastic can be bent with 3 prong pliers

Guaranteed*

Peace of mind for you and your patients

Optimum tooth contact

Superior retention

BPA free

Healthier patients

Gentle on tooth enamel


No scratching of the teeth

Stain resistant**

Plastic remains clear

Affordable

Contact us for pricing information



311 Viola St. | Oshawa, ON L1H 3A7

Phone: 905 436-3133

Toll free: 800 267-8463

Fax: 905 723-2331

Email: info@orthodont.ca

Website: www.orthodont.ca

* The plastic strap on the ClearBow assembly is guaranteed not to break or detach when the orthodontic appliance is used under normal conditions specified by the doctor. ** Beverages and food do not stain the ClearBow. It is important to properly clean the ClearBow on a regular basis; please visit clearbow.com for cleaning instructions.